**English**

**INSTRUCTIONS FOR USE**

For professional use only

**Español**

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Solo para uso profesional

**Français**

**MODE D'EMPLOI**

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Connexio Self-Etch Bond is a dual-curing, self-etching bonding agent reinforced with nano-particles for creating a durable bond between tooth substance and light-, self- and dual-curing restoratives without marginal leakage. This all-in-one adhesive consists of 2 components, which eliminates the need for a separate etching step, yet still achieves the same adhesive properties as total-etch bonds. Connexio Self-Etch Bond tolerates residual moisture and can be used with a variety of light, self and dual-curing composites (see contraindications).

**INDICATIONS:**

- Direct self- or dual-curing composite restorations and core buildups
- Direct light-curing composite-based restorations
- Fixing root pins/posts with dual- or self-curing composite cements
- Indirect restorations; with the application of dual-curing/self-curing composite cements for bonding inlays, onlays, crowns and bridges

**CONTRAINDICATIONS:**

- 1.** Connexio Self-Etch Bond contains organic acids, BIS-GMA, HEMA, TMPTMA, campherchinon, amines (DABE), BHT, catalysts, fluorides and ethanol. Do not use in patients with known allergies to these ingredients.
- 2.** Check instructions for use with any core buildup material for additional information.
- 3. Warnings and precautionary measures:**
  - Phenolic substances (e.g. temporary materials containing eugenol) disturb polymerization and therefore should not be used with Connexio Self-Etch Bond.
  - Avoid contact with the oral mucosa. Contact of the oral mucosa with Connexio may result in a whitish discoloration caused by the coagulation of proteins. This is a reversible occurrence that subsides within a few days. In case of contact with the eyes, rinse with plenty of clean water and consult a physician.
  - Application brushes provide precise and thin application of Connexio Self-Etch Bond.

**INSTRUCTIONS FOR USE:**

- 1. PREVENTING CONTAMINATION** – A dry working field is essential due to the subsequent application of the restorative. Avoid contamination with blood or saliva. Using a rubber dam is recommended.
- 2. PREPARATION** – Prepare the cavity according to the principles of the adhesive restoration technique. Rinse cavity thoroughly with water. Remove excess moisture with a faint air jet. Do not over dry dentin. The dentinal surface of the cavity should remain slightly moist. If the tooth has not been prepared (e.g. cervical defect) the surface must be thoroughly cleaned and enamel margin roughened with a finishing diamond. Areas in close proximity to the pulp should be protected with a suitable calcium hydroxide-based liner.
- 3. MIXING** – Connexio Self-Etch Bond is dual-curing and therefore must be applied immediately after mixing. Mix 1 drop of Part 1 and 1 drop of Part 2 on a mixing palette with an application brush for approximately 2 sec (produces a self-etching adhesive).
- 4. CONNEXIO SELF-ETCH BOND APPLICATION** – Apply the adhesive in a layer of medium thickness to the enamel/dentin and rub into the tooth surface for 20 sec.

**CAUTION:** The mixed bonding is light curing and thus an intensive exposure to ambient light should be avoided. The OP-Light during the application should be dimmed.

**5. DRYING** - Dry the adhesive layer for at least 5 sec with an air syringe.

**6. POLYMERIZATION –**

**a. Direct Restoration**

Light, self and dual-curing materials: The bonding layer must be polymerized with blue light (halogen or LED lights) for 10 sec when using light, self and dual-curing restoratives.

**b. Indirect restoration**

Light-curing luting composites: The bonding layer must be polymerized with blue light (halogen or LED lights) for 10 sec when using light-curing luting materials.

Self-curing or dual-curing luting materials: When using dual or self-curing luting materials, Connexio Self-Etch Bond also provides excellent adhesion without light-curing. Connexio Self-Etch Bond can also be light-cured for 10 sec with blue light (halogen or LED lights) as an optional step.

Luting root posts:

(a) Massage the mixed Connexio Self-Etch Bond (Part 1 and 2) into the occlusal surface of the core and root canal for 20 sec with an applicator. Disperse solvent with an oil-free air-jet for 5 sec. Remove excess material by dabbing it off with an absorbent strip of paper. Do not light-cure the bonding layer.

(b) Having prepared the root post in accordance with the manufacturer’s instructions, insert it into the root canal using a suitable luting material. The aim is to obtain a slight excess of cement.

(c) Light-cure for at least 40 sec to fix the root post and cure the Connexio Self-Etch Bond applied to the occlusal surface.

**7. ADVICE ON LUTING MARYLAND BRIDGES** – For luting Maryland bridges, the enamel must either be roughened or conditioned with a phosphoric acid etching gel before applying Connexio Self-Etch Bond.

Self-etching bonding materials are less effective on unprepared enamel. Applying bonding material to enamel that has not been ground or beveled can lead to margin discoloration. For optimum adhesion results with Connexio Self-Etch Bond, the enamel should be prepared in advance in class III, IV, V cavities and with diastema closures.

**Advice for storage and usage:**

Tightly close Connexio Self-Etch Bond bottles after use. Store in an upright position to guarantee flowback of the liquid. Store at temperatures between 39°F - 73°F (4°C - 23°C). Allow Connexio Self-Etch Bond to reach room temperature before use. Do not expose Connexio Self-Etch Bond to direct sunlight or the chair light.

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#### SAFETY DATA SHEET (SDS)

##### SECTION I - PRODUCT/COMPANY IDENTIFICATION

**REF 310270 CONNEXIO SELF-ETCH BOND, PART A AND B**

Application: Two-part, self-etching bonding liquid for dental bonding  
Company: CENTRIX, INC., 770 River Road, Shelton, CT 06484  
Toll-Free: 800-235-5862 (USA & Canada)

##### SECTION II - COMPOSITION/DATA ON COMPONENTS

Designation of Risk: Xi, irritant.    
Special Details of Risk (R-categories): 36-38  
Detrimental effects/symptoms: Irritant in case of eye, skin and inhalation contact. Harmful if swallowed. May cause sensitization by skin contact.  
**PART A: HAZARD CLASSIFICATION:** Corrosive  
NFPA ratings (scale 0 - 4) Health = 2 Fire = 0 Reactivity = 0  
HMIS-ratings (scale 0 - 4) Health = 2 Fire = 0 Reactivity = 0  
**PART A: HAZARD CLASSIFICATION:** Flammable  
NFPA ratings (scale 0 - 4) Health = 0 Fire = 3 Reactivity = 0  
HMIS-ratings (scale 0 - 4) Health = 0 Fire = 3 Reactivity = 0

##### SECTION III - COMPOSITION/DATA ON COMPONENTS

**PART B** Description: Mixture of the substances listed below with nonhazardous additions.

<b>Components:</b>	%	<b>CAS NO.</b>
Acidic adhesive monomer	50-100%	1565-94-2
BIS-GMA	5-10%	868-77-92
hydroxyethyl methacrylate	5-10%	112813-65-7

##### PART B

Description: Mixture of the substances listed below with nonhazardous additions.

<b>Components:</b>	%	<b>CAS NO.</b>
ethanol	50-100%	64-17-5
Initiator	2.5-5%	- - -

Exact % content withheld due to trade secret restrictions.

##### SECTION IV - FIRST AID MEASURES (APPLIES TO BOTH PART A AND B)

Eye Contact:	Flush immediately with plenty of water. Consult a physician if irritation persists
Skin Contact:	Wash with plenty of water and soap.
Ingestion:	Consult a physician.
Inhalation:	Supply fresh air; consult doctor in case of complaints.

##### SECTION V - FIRE-FIGHTING MEASURES (APPLIES TO BOTH PART A AND B)

Suitable Extinguishing Media: Foam, dry powder, carbon dioxide, water spray  
Unusual fire and explosion hazards: None  
Specific intervention methods: Appropriate protective equipment must be worn in case of fire.  
**PART B ONLY:** For safety reasons unsuitable extinguishing agents: Water with full jet

##### SECTION VI - ACCIDENTAL RELEASE MEASURES (APPLIES TO BOTH PART A AND B)

Person-related safety precautions: Wear protective equipment. Keep unprotected persons away.  
Measures for environmental protection: Do not allow to enter sewers/ surface or ground water.  
Dilute with plenty of water.  
Measures for cleaning/collecting:  
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust). Use neutralizing agent. **PART A ONLY:** Ensure adequate ventilation.

##### SECTION VII - HANDLING AND STORAGE (PART A AND B)

###### Handling

Information for safe handling: For dental use only. Protect against electrostatic charges.  
**PART B ONLY:** Keep ignition sources away - Do not smoke. No special measures required.

###### Storage

Keep receptacle tightly sealed.

##### SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

###### PART A:

Additional information about design of technical systems: No further data; see item 7.  
Components with limit values that require monitoring at the workplace:  
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.  
Personal protective equipment: General protective and hygienic measures:  
Keep away from foodstuffs, beverages and feed.  
Wash hands before breaks and at the end of work.  
Avoid contact with the eyes and skin.  
Breathing equipment: Not required.  
Protection of hands: Protective gloves  
Eye protection: Safety glasses

##### SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

###### PART A:

Form: Fluid	Color: Colorless	Odor: Characteristic
Flash point: Not applicable.		
Density at 20°C (68°F):	1.2 g/cm³	
Solubility in / Miscibility with Water: Partly miscible.		
pH-value at 20°C (68°F):	1.5	

###### PART B:

Form: Fluid	Color: Colorless	Odor: Characteristic
Change in condition		
Flash point: 21°C (70°F)		
Auto igniting: Product is not selfigniting.		
Danger of explosion: Product is not explosive. However, formation of explosive air/vapor mixtures are possible.		
Density at 20°C (68°F): 0.87 g/cm³		
Solubility in / Miscibility with Water: Fully miscible.		
Components with limit values that require monitoring at the workplace:		
64-17-5 ethanol		
PEL 1900 mg/m³, 1000 ppm		
REL 1900 mg/m³, 1000 ppm		
TLV Short-term value: 1880 mg/m³, 1000 ppm		
Other Data: Not available or none.		

##### SECTION X - STABILITY AND REACTIVITY (APPLIES TO BOTH PART A AND B)

No dangerous reactions or decomposition products known.

##### SECTION XI - TOXICOLOGICAL INFORMATION

###### PART A:

Acute toxicity:  
LD/LC50 values that are relevant for classification: Acidic adhesive monomer  
Oral LD50 > 2000 mg/kg (rat)  
Primary irritant effect on the skin: Irritant to skin and mucous membranes.  
Primary irritant effect on the eye: Irritating effect.  
Additional toxicological information: The product shows the following dangers according to internally approved calculation methods for preparations: Corrosive  
**PART B:**  
Acute toxicity:  
Primary irritant effect on the skin or eyes: No irritant effect.  
Sensitization: No sensitizing effects known.

##### SECTION XII - ECOLOGICAL INFORMATION (APPLIES TO BOTH LIQUIDS A AND B)

Water hazard class 1 (Self-assessment): slightly hazardous for water. Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

##### SECTION XIII - DISPOSAL CONSIDERATIONS (APPLIES TO BOTH LIQUIDS A AND B)

Must not be disposed of together with household garbage. Do not allow product to reach sewage system. Dispose in accordance with federal, state and local regulations

##### SECTION XIV - TRANSPORT INFORMATION

**PART A:** DOT regulations: Label 8

###### PART B:

<b>DOT regulations</b>		
Hazard class: 3	Identification number: UN1170	
Packing group: II	Proper shipping name (technical name): ETHANOL, mixture	
Label 3		
<b>Land transport ADR/RID (cross-border):</b>		
ADR/RID class: 3 Flammable liquids		Danger code (Kemler): 33
UN-Number: 1170		Packaging group: II
Label: 3		
Description of goods: 1170 ETHANOL (ETHYL ALCOHOL), mixture		
<b>Maritime transport IMDG:</b>		
IMDG Class: 3	UN Number: 1170	Label 3
Packaging group: II	EMS Number: F-E,S-D	Marine pollutant: No
Proper shipping name: ETHANOL (ETHYL ALCOHOL), mixture		
<b>Air transport ICAO-TI and IATA-DGR:</b>		
ICAO/IATA Class: 3	UN/ID Number: 1170	
Label 3	Packaging group: II	
Propper shipping name: ETHANOL, mixture		

##### SECTION XV - REGULATORY INFORMATION

Sara: Sections 315, 355; none of the ingredients are listed.  
TSCA: all of the ingredients are listed  
Proposition 65: none of the ingredients are listed  
Canadian Ingredient Disclosure list (limit 1%) 868-77-9 2-hydroxyethyl methacrylate (PART A)  
Product related hazard information: The product has been classified and marked in accordance with directives on hazardous materials.  
Hazard symbols PART A: Corrosive  
Hazard symbols PART B: Flammable  
Hazard-determining components of labelling: 2-hydroxyethyl methacrylate, Acidic adhesive monomer (PART A)  
Risk phrases: Causes burns (PART A) Risk phrases: Flammable (PART B)  
May cause sensitization by skin contact (PART A)  
Safety phrases: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water. Wear suitable gloves and eye/face protection.  
Keep away from sources of ignition - No smoking (PART B)

##### SECTION XVI - OTHER INFORMATION

Product handling exclusively by dental specialists per instructions for use.

##### PRODUCT FOR PROFESSIONAL USE ONLY.

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship. Department issuing MSDS: Scientific department  
Contact: Tel.: +49 (0)4721-39 99 6-0

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© 2016 Centrix Inc. 012091A REV E Date Revised: 07/2016