# couture<sup>TM</sup> Diamond Polishers


couture<sup>TM</sup>  
Advanced Finishing & Polishing

## ENGLISH • INSTRUCTIONS FOR USE

### Couture™ Diamond Polishers

REF 390001 – Diamond Finishing and Polishing Kit  
REF 390002 – Diamond Finishing Cups Medium  
REF 390003 – Diamond Finishing Points Medium  
REF 390004 – Diamond Polishing Cups UltraFine  
REF 390005 – Diamond Polishing Points UltraFine

#### PRODUCT DESCRIPTION

- Unique polishing cups and points with highly engineered core designs so that the Polishers retain their shapes during polishing procedures.
- UltraFine (blue) and Medium (green) Polishers are highly loaded with ultrafine diamonds for superior finished surfaces providing more life-like reflectance and luster.
- For use in right angle, latch style handpieces at low speeds (8,000-20,000 rpm).
- Length of the working section: Point: 9 mm; Cup: 7 mm
- Maximum diameter of the working section: Point: 3 mm; Cup: 6 mm
- Autoclavable for multiple uses

#### INDICATIONS FOR USE

- For initial finishing and final high gloss polishing of convex and pronounced concave surfaces of esthetic/composite restorations
- Maintenance of existing esthetic/composite restorations

#### PRACTICAL TIPS

- Use Couture Polishers on composite restorations after pre-polishing with rotary diamond instruments (15 µm yellow).
- Medium (green) Polishers enable easy shaping and well prepared surfaces.
- UltraFine (blue) Polishers enable final, ultra-high gloss surfaces, bringing out the maximum esthetics of any composite material.
- Use copious water spray during the entire procedure (minimum 50 mL/min).
- Apply light force to instruments, not exceeding 0.5 lbs (2N) of pressure.
- Use slow speeds between 8,000 and 20,000 rpm maximum (note: higher speeds will cause more rapid deterioration of Polisher and rougher surfaces).

#### MAINTENANCE AND STERILIZATION

- Instruments are delivered non-sterile. We recommend sterilization prior to first use and then cleaned, disinfected and sterilized after each use.

**Pre-cleaning:** Remove all dirt on surfaces using a paper towel.

#### Automatic cleaners/disinfectors:

- 1) Place instruments in the machine.
- 2) Set the machine on standard cycle and clean for at least 5 minutes.
- 3) Rinse for at least 3 minutes.

**Disinfection:** Disinfection solutions may be used, follow label instructions. For automatic cleaners/disinfectors, 10 minutes rinsing at 175°F (80°C) may be performed as a last step.

**Sterilization:** For sterilization processes, we refer to the ISO 17664 standard and recommend the following procedures:

- 1) 273°F/134°C Cycles:  
  - Tmin. = 134°C Tmax. = 138°C
  - Pressure = 45 PSI/3.15 bar (abs)
  - Time = 4 minutes
- 2) 250°F/121°C Cycles:  
  - Tmin. = 121°C Tmax. = 125°C
  - Pressure = 30.5 PSI/2.10 bar (abs)
  - Time = 16 minutes

**Maintenance:** Discard any blunt or damaged Polisher. Regularly check instruments to ensure that instruments are not damaged or worn; that edges and corners are not notched and remain regular in shape.

#### RISK WARNINGS

- Insufficient water cooling may lead to irreversible damage to restoration and adjacent tissues.
- Temperatures above 280°F (140°C) should be avoided as they could increase the risk of mandrel fracture during usage.
- Avoid heavy operating force exceeding 0.5 lbs (2N) of pressure due to increased risk of fracture of the plastic shaft and of the tip of the Point Polisher.
- To ensure traceability of the instruments, we recommend keeping the original packaging.

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