ENGLISH INSTRUCTIONS FOR USE

Couture™ Diamond Polishing Brushes

REF 390001 - Diamond Finishing and Polishing Kit

REF 390006 - Diamond Polishing Brushes

PRODUCT DESCRIPTION

- Unique cellulose brush with filaments hyper-filled with ultra-fine diamond particles for high gloss polishing of esthetic/composite restorations.
- Only short polishing times are needed to obtain optimal results.
- Unique cylindrical shape, with more than 600 extra-long filaments.
- When used dry, filaments retain firmness for more aggressive polishing.
- Filaments change from hard to ultra soft with water absorption. As it absorbs water, the softened filaments release the ultra-fine diamond particles.
- The free ultra-fine diamond particles adapt to the tooth surface structure and significantly increase polishing performance compared to currently available results with other brushes.
- When wet and softened, the brush has a low-vibration and is atraumatic to the soft tissues, making it safe for even lapping onto gingivae.
- Stable plastic mandrel for use in latch-type, contra-angle handpieces.
- Autoclave instrument before first use.
- Single use (due to water and debris absorption, possible deformation).

INDICATIONS FOR USE

- Final high-gloss polishing of all esthetic/composite restorations
- Maintenance of existing restorations
- For concave and convex tooth surface morphologies

PRACTICAL TIPS

- Autoclave instruments before first use.
- Use Couture Diamond Polishing Brush after pre-polishing with rotary diamond instruments (15µm yellow).
- Use Couture Diamond Polishing Brush on the dental surface starting with hard filaments for 15 seconds under light pressure (0.5 LBS/2 N) and low speed (5000 RPM).
- Subsequently add water spray for few seconds.
- The stiffness of the filaments gradually change to ultra soft. Adjust spray if there is too much water.
- Continue polishing at low speed for about 15 seconds more without pressure, with only light wiping movements to achieve final high gloss.
- Apply Couture Diamond Polishing Brush without paste.
- Single use (due to water and debris absorption, deformation).
- Couture Diamond Polishing Brush can be safely disposed as it is made entirely of recyclable materials.

STERILIZATION

- Sterilization should be done according to validated methods. If possible, use pre- or post-vacuum autoclave Class B with vacuum drying.
- For sterilization processes, we refer to the ISO17664 standard and recommend the following procedures:

- 1) 273°F/134°C Cycles:
 - Tmin. = 134°C
 - Tmax. = 138°C
 - Pressure = 45PSI/3.15 bar (abs)
 - Time = 4 minutes

- 2) 250°F/121°C Cycles:
 - Tmin. = 121°C
 - Pressure = 30.5 PSI/2.10 bar (abs)
 - Time = 16 minutes

RISK WARNINGS

- Single use (due to water and debris absorption and deformation).
- Additional sterilization cycles may also change the plastic structure and increase the risk of fracture of the shaft. Check after each sterilization cycle.
- Temperatures above 280°F (140°C) should be avoided as the plastic shaft becomes increasingly brittle and subject to fracture during subsequent use.
- Avoid contact pressure of over 0.5 LBS/2N, because of fracture risk of the plastic shaft.
- To ensure traceability of the instruments, we recommend keeping the packaging.

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Couture™ Diamond Polishing Brushes

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