BUILD-UP RIPRIMITIVO ENCORE®D/C MINIMIX®

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ENCORE® D/C MINIMIX® CORE BUILDUP
Dual-cure composite resin
Indicated for use as a core buildup composite under crowns and bridges. May also be used to cement posts when used in conjunction with a suitable dentin adhesive.
INSTRUCTIONS FOR USE
1. Isolate the tooth to be prepared. We recommend the use of a rubber dam or cotton rolls.
2. Remove all existing restorations and decay from the tooth. At this time you may place any pins or posts if necessary.
3. Acid etch the entire preparation for 20 seconds using (32% to 40%) phosphoric acid. We recommend the use of Centrix Onyx™ 40% phosphoric etch gel.
4. Rinse with a water spray for 30 seconds and lightly air dry. The enamel of the preparation should have a frosted appearance.
5. Place bonding system of your choice onto the preparation, according to the manufacturer’s instructions. Do not use one bottle bonding systems. When placing the primer, be sure that all of the solvent has been evaporated and the preparation has a glossy appearance. Before light curing, the bonding agent should be thinned with a light air spray.
6. Build up core

- Insert Encore D/C Mini-Mix cartridge into a Centrix Snap-Fit™ Syringe.
- Remove cap from MiniMix cartridge and attach mixing nozzle.
- Purge a small amount from the tip to assure proper mixing and flow.
- Inject the material into the prep, keeping the needle at the bottom of the prep to reduce the chance of voids, fill from the bottom up.

IF USING A MATRIX BAND

- To facilitate placement, a matrix band may be placed around the prepared tooth.
- Place the mixing tip directly into the preparation and squeeze the trigger of the applicator gun to express the paste.
- If a matrix band was used, test with an explorer to insure that material has set prior to removal. Once set, cut and finish the preparation, using water-cooled diamonds or burs.

IF USING A PREFABRICATED CROWN FORM

- Place the mixing tip directly into the crown form and squeeze the trigger of the cartridge to express the paste.
- Place the loaded crown form onto the prepared tooth and hold into place until the paste is fully set.

WORKING TIME				
	Working Time	Self-CureTime	Light-Cure	Total Setting Time
Self-Cure	1:30*	3:30*		5 minutes
Light-Cure	1:30* or less		:40"	2:10* or less

Encore D/C MiniMix core buildup paste has a 1:30 minute working time. The material may be manipulated using a composite instrument. A flat-ended interproximal carver (IPC) is recommended. Light activate with a standard dental halogen lamp or LED curing unit for 40 seconds (per 4mm layer). A thinner layer (<2mm) will need only a 20 second light exposure. To complete the chemical cure activation, wait 5 minutes before cutting and finishing the prep. If a matrix band was used, test with an explorer to insure that material has set prior to its removal. Once set, cut and finish the preparation, using water-cooled diamonds or burs.

PRECAUTIONS

- If pulp is exposed, use a suitable liner.
- Contact with skin and eyes should be avoided. If accidental contact occurs, rinse immediately with plenty of water and consult a physician if necessary.
- Do not use with materials containing eugenol. This could inhibit polymerization.
- Water and oil-contaminated air must be avoided, as they inhibit polymerization at the contact area.
- Avoid contact with gingival tissue and adjacent teeth.
- Placement of any composite is contraindicated if a dry field is not possible.
- Do not use Encore D/C MiniMix core buildup paste for direct application on dentin near the pulp or if the pulp chamber has been exposed.
- There are no known systemic side effects to date, however, do not use in patients with known sensitivities to Bis-GMA, benzoyl peroxide or composite resin materials.

Store away from sunlight, at room temperature 23°C (72°F)

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SAFETY DATA SHEET (SDS)

In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision
SECTION I—IDENTIFICATION
Product Trade Name: ENCORE D/C MINIMIX Part (Item) Number: REF 310254, 310255, 310129, 310130
Product Name: Dental Core Build-Up material
Manufacturer: CENTRIX, INC. Address: 770 River Road, Shelton, CT 06484 U.S.A.
Emergency or Information Telephone Numbers: (203)929-5582, Toll-Free: (800)235-5862 (USA & Canada) Fax: (203)929-6804
Date Revised: 3/2015

SECTION II—HAZARDS IDENTIFICATION
Description of Risk: Xi, Irritant
Special Details of Risk (R-categories): 36-38 Detrimental effects/ symptoms; Irritant in case of eye, skin and inhalation contact.
Harmful if swallowed. May cause sensitization by skin contact.

SECTION III—PHYSICAL DATA
Hazardous Components:
Blend of multifunctional methacrylates:

%	CAS NO.	SYMBOL	R	OTHER LIMITS
26-32	N/A	Xn/Xi	22/36/37/38	

Initiator:

%	CAS NO.	SYMBOL	R	OTHER LIMITS
0.10 - 0.35	94-36-0	E/Xn	2/7/36/43	

Exact % content withheld due to trade secret restrictions.

SECTION IV—FIRST AID MEASURES
Risk: Immediately wash skin with soap and water. If irritation occurs seek medical advice.
Ingestion: If swallowed do not induce vomiting. Give large quantities of water and milk. Seek medical attention.
Inhalation: Remove to fresh air. If breathing becomes difficult, contact a physician.
Eyes: Flush with plenty of water. Seek immediate medical attention by an ophthalmologist.

SECTION V—FIRE-FIGHTING MEASURES.
The product is not considered to be a fire hazard. The information below is given for extinguishing a fire generated by other sources.
Suitable extinguishing media: Water, Chemical Foam, Carbon Dioxide, or Dry Chemical.
Forbidden extinguishing media: Not determined.
Special firefighting measures: None. Wear self-contained breathing apparatus.
Unusual fire and explosion hazards: Heat can cause polymerization with rapid release of energy.
Special protection equipment: Self Contained Breathing Apparatus if toxic fumes are generated.

SECTION VI— ACCIDENTAL RELEASE MEASURES
Personal precautions: Adopt some precautions of section VII.
Environmental precautions: Keep spilled material out of sewers.
Reclaiming methods: Absorb spills with inert material – keep spilled material out of sewers.

SECTION VII— HANDLING AND STORAGE
Steps to be taken in case material is released or spilled: Absorb with inert materials such as dry sand and place in closed container for disposal as solid waste.
Waste Disposal Method: This disposal method is recommendation only. Be sure to dispose of materials in accordance with all local, state and federal regulations applicable.

SECTION VIII— EXPOSURE CONTROLS / PERSONAL PROTECTION
Respiratory Protection, Ventilation, Protective Gloves, Eye Protection, and other protective equipment:
Not required under normal condition of use.

SECTION IX— PHYSICAL AND CHEMICAL PROPERTIES
Appearance: Colored/uncolored paste
Explosion limits: Not determined
Odor: None
Water solubility: Insoluble
Boiling point: Not determined
Evaporation rate: Not determined
Flash point: Not determined
Other Data: Not available or none.

SECTION X—STABILITY AND REACTIVITY
Stability: Stable if stored as directed.
Conditions to avoid: Heat, light, aging and contamination.
Materials to avoid: Reducing and oxidizing agents, peroxides and amines.
Hazardous decomposition products: Oxides of Carbon.
Hazardous polymerization products: May occur.

SECTION XI—TOXICOLOGICAL INFORMATION
Ingestion: Uncured material may be harmful or fatal if swallowed.
Eye contact: May cause irritation and damage if not removed promptly.
Skin contact: Prolonged or repeated exposure to uncured material may cause irritation or skin rash.
Inhalation: Prolonged or excessive inhalation may cause respiratory tract irritation.
Prolonged exposure: Not applicable.

SECTION XII— ECOLOGICAL INFORMATION
There is no specific information on the product. Use according to good working practice.
Avoid releasing the product into environment.

SECTION XIII— DISPOSAL CONSIDERATIONS
Unpolymerized (uncured) material may be hazardous waste.
Incinerate uncured material and dispose in accordance with local regulations.

SECTION XIV— TRANSPORT INFORMATION
Additional Information: Not dangerous cargo.

SECTION XV— REGULATORY INFORMATION
CLASSIFICATION IN ACCORDANCE WITH EEC DIRECTIVES.
Hazard Classification Symbol: Xn
R22 Harmful if swallowed.
R36-38 Irritating to eyes, respiratory system and skin.
S36 Wear suitable protective clothing.
Ingredients determining the label: Blend of multifunctional methacrylate

SECTION XVI— OTHER INFORMATION
Product handling is exclusively by dental specialists per instructions for use.

CAUTION
The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user’s responsibility.

Encore® D/C