10°C -24°C / 50°F -75°F

ENGLISH

INSTRUCTIONS FOR USE

Intermediate Restorative Material

30 x .20 pre-filled syringe tips

REF 310064 with fluoride **REF 310065** without fluoride

Contra-indications: **Do not use in presence of eugenol**

- 1) Clean and disinfect prep. On exposed pulp, place a calcium hydroxide pulp-capping agent (Calciject™). On vital dentin, a desensitizing liner (Glu/Sense™) is recommended. Do not acid etch or use bonding agents.
- 2) Prep should be as dry as possible.
- 3) Place TEMPIT® ULTRA tip in syringe. Remove cap and inject. Use a matrix when indicated. Shape and contour with conventional composite instruments. To reduce stickiness, dip instruments in alcohol or refrigerate Tempit® Ultra tips before use.
- 4) Light activate with a conventional halogen light source for 40 seconds. TEMPIT® ULTRA will cure up to 6 mm deep.
- 5) To adjust occlusion or remove TEMPIT® ULTRA after it is set, use water-cooled diamonds or burs. Keep away from sources of strong visible light. Do not freeze.

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SAFETY DATA SHEET (SDS)

In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

SECTION I - PRODUCT/ COMPANY IDENTIFICATION

Product Trade Name:
Tempit® Ultra-F - Light-Cure Temporary Restorative Paste with Fluoride
Tempit® Ultra - Light-Cure Temporary Restorative Paste
Synonyms: Light-Cure Temporary Restorative with/without Fluoride Paste
Recommended Use: Dental Material
Part (Item) Number: 310064, 310065
Manufacturer: CENTRIX INC, 770 RIVER ROAD, SHELTON, CT 06484 U.S.A.
Emergency Telephone Number: (800) 235 - 5862
Information Telephone Number: (203) 929-5582

SECTION II – HAZARDS IDENTIFICATION

For potential health hazards: SEE Section 11 below.
No other significant hazards to man or the environment are reasonably foreseen.

SECTION III - COMPOSITION / INFORMATION ON INGREDIENTS

Product Contents:
Tempit® Ultra-F - Light-Cure Temporary Restorative Paste with Fluoride
Tempit® Ultra - Light-Cure Temporary Restorative Paste
Hazardous Ingredients:
Ingredient HAZARDOUS DATA Hazardous Data Legend
Sodium Fluoride Solid Health: 3 4 - Extreme
CAS: 7681-49-4 Flammability: 0 3 - High
UN Designation: UN1690 Reactivity: 1 2 - Moderate
Class: 6.1. Toxic 1 - Slight
Packing Group: III 0 - Insignificant
Exact % content withheld due to trade secret restrictions.

SECTION IV - FIRST AID MEASURES

SKIN: Wash with soap and water.
INGESTION: Induce vomiting if large amount of material is swallowed.
INHALATION: No harm anticipated.
EYES: Flush with copious amounts of water.

SECTION V - FIRE-FIGHTING MEASURES

SUITABLE EXTINGUISHING MEDIA: Water contact acceptable.
FORBIDDEN EXTINGUISHING MEDIA: None.
UNUSUAL FIRE AND EXPLOSION HAZARDS: None specific.
SPECIAL PROTECTION EQUIPMENT: None.

SECTION VI - Accidental Release Measures

PERSONAL PRECAUTIONS: Avoid contact of skin or eyes. Wash hands with soap and water if contact occurs.
ENVIRONMENTAL PRECAUTIONS: Dispose of in accordance to all federal, state, local regulations.
RECLAIMING METHODS: Wipe up with soap, water and paper towels.

SECTION VII - HANDLING AND STORAGE

HANDLING: None. STORAGE: Store at temperatures not exceeding 75°F (24°C).
CONTAINERS: Avoid exposure to moisture and light.
INCOMPATIBILITIES: None specific. OTHER PRECAUTIONS: None.

SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

GENERAL MEASURES: None.
SKIN PROTECTION: Avoid prolonged contact skin; wash with soap and water.
EYE PROTECTION: Goggles/Eyeshields.
HANDS PROTECTION: Avoid prolonged contact skin; wash with soap and water.
LIMITS: Respiratory protection: Not required
OTHER EQUIPMENT:

SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

APPEARANCE: Tooth shade paste EXPLOSION LIMITS: N/A
ODOR: Negligible odor WATER SOLUBILITY: Virtually insoluble
BOILING POINT: N/A EVAPORATION RATE: Not established FLASH POINT: N/A

SECTION X - STABILITY AND REACTIVITY

STABILITY: Generally stable.
CONDITIONS TO AVOID: Light exposure, heat.
MATERIALS TO AVOID: N/A
HAZARDOUS DECOMPOSITION PRODUCTS: N/A - unless exposed to high temperatures and/or light and only when large quantities of the material are involved.

SECTION XI - TOXICOLOGICAL INFORMATION

POTENTIAL ACUTE HEALTH EFFECTS: Possibility of skin irritation or sensitization if not washed off promptly
POTENTIAL CHRONIC HEALTH EFFECTS: Possibility of skin sensitization or irritation upon prolonged or repeated exposure.
SIGNS AND SYMPTOMS OF EXPOSURE: Skin/eye irritation.
MEDICAL CONDITIONS GENERALLY AGGRAVATE: None known.
CARCINOGENIC OR POTENTIAL CARCINOGENIC: No ingredient listed as a carcinogen.
OTHER POTENTIAL HEALTH EFFECTS: None known.

SECTION XII - ECOLOGICAL INFORMATION

No ecological damage or impact anticipated when handled as indicated.

SECTION XIII -DISPOSAL CONSIDERATIONS

WASTE DISPOSAL: Do not dispose into sink. Dispose in accordance with all federal, state, and local regulations.

SECTION XIV - Transport Information

No hazardous goods within the meaning of transport regulations.

SECTION XV - REGULATORY INFORMATION

U.S. FDA DENTAL DEVICE: Class II product
EC MEDICAL DEVICE DIRECTIVE: Class IIa product

SECTION XVI - OTHER INFORMATION

Product handling exclusively by dental specialists per instructions for use.

CAUTION

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility

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EMERGO AUSTRALIA

Level 20, Tower II, Darling Park, 201 Sussex Street
Sydney, NSW 2000, Australia