



SAFETY DATA SHEET (SDS)

In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

SECTION I - PRODUCT / COMPANY IDENTIFICATION

Product Name: HeatSleeve™ Composite Warmer
Part Number: 330553
Description: Self-heating Dental Composite warmer
Company: Centrix, Inc.
Address: 770 River Road City, State, Zip: Shelton, CT 06484
Emergency / Information Telephone Numbers: 203-929-5582
Toll-Free: 800-235-5862 (USA & Canada) Date Revised: 11/2015

SECTION II – HAZARD IDENTIFICATION

Designation of Risk: Xi, irritant  
Special Details of Risk (R-categories): 36-38, 41-43
Detrimental effects/symptoms: Irritant in case of eye, skin and inhalation contact. Risk of damage to eyes. Harmful if swallowed. May cause sensitization by skin contact.

SECTION III - COMPOSITION / INFORMATION ON INGREDIENTS

Ingredient	CAS-No.	% by weight
Bronopol	52-51-7	<0.05
5-chloro-2-methyl-2H-isothiazol-3-one	26172-55-4	<0.005
Sodium Acetate	127-09-3	85-95
Distilled Water	7732-18-5	5-15
Dipropylene Glycol	25265-71-8	0.01
Methylisothiazolinone	2682-20-4	<0.005

Exact % content withheld due to trade secret restrictions.

SECTION IV - FIRST AID MEASURES

Inhalation: If cough, shortness of breath or other breathing problems occur, move to fresh air
Skin Contact: Wash promptly with copious amounts of water.
Eye Contact: Wash promptly with copious amounts of water.
Ingestion: Consult a physician.

SECTION V - FIRE-FIGHTING MEASURES

Suitable Extinguishing Media: CO2, powder or water spray
Extinguishing Media which should not be used: Not applicable
Special exposure hazards from combustion products: Not applicable
Special protective equipment for fire fighters: Not applicable

SECTION VI - ACCIDENTAL RELEASE MEASURES

Personal Precautions: Refer to section 8 of SDS for personal protection details.
Turn leaking containers leak-side up to prevent the escape of liquid.
Environmental precautions: Do not discharge into drains or rivers. Contain the spillage using bunding.
Clean-up procedures: Absorb into dry earth or sand. Transfer to a closable, labelled salvage container for disposal by an appropriate method.

SECTION VII - HANDLING AND STORAGE

Handling Precautions: Avoid contact with skin, eyes, and clothing; ensure good ventilation.
Safe Storage Conditions: No special requirements
OTHER RESTRICTIONS OR RECOMMENDATIONS: No further relevant information available

SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

Respiratory Protection & Ventilation: None required.
Protective Clothing or Equipment: Protective gloves
Unusual Work/Hygienic Practices: No special requirements

SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Colored fluid Odor: Odorless pH: not determined
Boiling Point/Range: undetermined Melting Point/Range: undetermined
Flash Point: Not applicable Auto-Ignition Temperature: Not applicable
Flammability: Not applicable Explosive Limits in Air: Not Determined
Solubility in Water: Soluble Reactivity in Water: Fully Miscible Other Data: None

SECTION X - STABILITY AND REACTIVITY.

Stability: Stable
Conditions of Instability: Stable under normal conditions.
Incompatibility with Various Substances: No relevant information available
Hazardous Decomposition Products: None known.
Hazardous Polymerization: Will not occur.
Conditions to Avoid: Not applicable

SECTION XI - TOXICOLOGICAL INFORMATION.

Potential Acute Health Effects: Damage to skin and eyes.
Potential Chronic Health Effects: None known.
Signs and Symptoms of Exposure: Skin and eye irritation, skin discoloration.
Medical Conditions Generally Aggravated by Exposure: None known.
Carcinogenic or Potential Carcinogenic Effects: None anticipated.
Other Potential Health Effects: None known.

SECTION XII - ECOLOGICAL INFORMATION.

None known to be applicable to the subject product or similar materials.

SECTION XIII - DISPOSAL CONSIDERATIONS.

Waste Disposal: Dispose in accordance with the local, state, and federal regulations.

SECTION XIV - TRANSPORT INFORMATION.

Additional Information: Not dangerous cargo.

SECTION XV - REGULATORY INFORMATION

Dental Device: Class I medical device for FDA and EC Medical Device Directive.

SECTION XVI - OTHER INFORMATION

Product handling is exclusively by dental specialists per instructions for use.

CAUTION: The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.

centrix®

Making Dentistry Easier.™

CENTRIX INCORPORATED
770 River Road, Shelton, CT 06484 USA
(203) 929-5582 (800) 235-5862 FAX (203) 929-6804
www.centrixdental.com



EC | REF | Dental-Contact Vertriebs OHG, Im Spiet 101, DE 26506 Norden

EMERGO AUSTRALIA
Level 20, Tower II, Darling Park, 201 Sussex Street
Sydney, NSW 2000, Australia

Instructions for Use


1



2



3



4

