



## BIOACTIVE. MOISTURE FRIENDLY. MARGIN INTEGRITY. PROVEN FLUORIDE RELEASE.



### FLUORIDE RELEASE

Champ® Pit & Fissure Sealant effectively releases fluoride under aqueous conditions over several weeks. As adjacent enamel is challenged, this fluoride may be available to help protect the teeth from decay. Our experience in fluoride release testing suggests that the fluoride release shown here would provide clinically beneficial attributes.

CHAMP PIT & FISSURE SEALANT CUMULATIVE FLUORIDE RELEASE DATA<sup>1</sup>



<sup>1</sup>Internal study, summary available upon request