


centrix[®]
MAKING DENTISTRY EASIER[®]

PERFECT
PREVENTION

EFFICIENT
DELIVERY

ACCURATE
IMPRESSIONS

**BEAUTIFUL
RESULTS!**



Celebrating 50 Years

"Like most practicing dentists, I experienced problems with products that made the delivery of quality dentistry more difficult than it needed to be. My first products, AccuDose® tubes and the C-R® syringes, changed the use of dental materials and made dentistry easier. 50 years later, I am proud that Centrix continues to introduce products and ideas that fulfill our mission."

Making Dentistry Easier.

~ Dr. William Dragan



MAKING DENTISTRY EASIER!

Today, Centrix manufactures and engineers hundreds of dental products focused primarily on general restorative and preventive dentistry, the heart and soul of most dental practices. As a leader in materials and delivery systems, Centrix continues to introduce inventive dental systems certain to become the new standard in dentistry.

OUR FIRST INNOVATIONS!



Mark IIIp™ Gun & Snap-Fit™ Gun

Easy, direct-to-prep delivery of today's dental materials. Made in the USA and Used in every dental office!



Benda® Brush Applicators

One-piece, bendable head, precision brushes. Made in the USA. Can you imagine dentistry without them?



AccuDose® & NeedleTube®

Direct-to-prep done better! Dispense exactly the material needed with no waste or cross-contamination. Centrix Tubes & Plugs are made in the USA and truly changed the industry!

OUR NEWEST INNOVATIONS!



NoCord™ VPS

Industry-First: Hemostasis and retraction built-in! Award-winning, one-step, self-retracting, impression material. Better marginal detail and better fluids control.



Encore® D/CB + Connexio®

Making post-endo treatment easier and more assured. For restorations that last! Increased bond strengths reduce polymerization shrinkage, virtually eliminating micro-leakage.



PREVENTION FOR LIFE®

Caries rates in kids are decreasing! However, caries rates in adults are increasing, rapidly, as patients decline to pay for treatment not covered by insurance – or are not offered preventive care, as a matter of course.

The Centrix Prevention for Life program provides tools and information to help offices update their preventive treatment protocols and how they discuss them with patients. **Webinars and in-office Lunch & Learn** programs that provide CE credits. **Risk Assessment Tools** that help you quickly, easily and visually share with your patients your diagnosis. **Prescription and Treatment Program** that increases patient health (as well as practice health). **Products** that are easy to place and your patients will love!

At Centrix, we believe that structured Risk Assessment quickly results in diagnosis of risk and a prescription for that patient (rinse, varnish, etc.), treatment plan and frequency for each patient, and easy patient acceptance!

Start a Prevention For Life® program today! Visit cariesprevention.com or call 800.235.5862



ADA Center for Evidence-Based Dentistry (November 2013) Topical Fluoride for Caries Prevention. Retrieved from <http://www.ada.org>

ORAL CARE AND PREVENTION

Centrix makes prevention easier for you and your patients



FluoroDose®

Award-Winning Fluoride Varnish with Xylitol



- Tastes Better! Feels Better! Better Priced!
- Application takes less than a minute.
- Unit-dose delivery means easier use, easier cleanup.

TOP FLUORIDE VARNISH, 6 YEARS IN A ROW!



DENTAL ADVISOR TOP AWARD: FLUORIDE VARNISH



LolloCaine®

Topical Anesthetic | 20% Benzocaine Gel

- 4 Great flavors: Bubble Gum, Cherry, Mint, Piña Colada.
- Rapid onset with no alcohol aftertaste!
- Unit-dose LolloPack® saves time and reduces cross-contamination.
- The LolloPack is the perfect “take-home pack” for post-op discomfort.



E-Lieve™

Vitamin E Treatment | Made in the USA

E-Lieve helps to relieve:

- Chapped lips before, during and after dental procedures.
- Dryness and irritation when using rubber dams.
- Sore tissue from deep scaling procedures.
- Discomfort from wearing new dentures.
- Irritation caused by orthodontic brackets.

A Faster, Easier Solution for Moisturizing and Soothing Minor Soft Tissue irritation.



Champ®

Pit & Fissure Sealant | Glass-Filled Acrylate Resin

- **Bioactive:** hydrophilic resin forms a strong resin-hydroxyapatite complex that seals out micro-leakage.
- **Moisture-friendly formula:** tolerates moisture making it easier to bond to the tooth.
- **Margin Integrity:** lightly filled (37%) resin is durable and long-lasting; adapts to and seals the margins, helping to eliminate micro-leakage.



Exposé®

Caries Indicator on a Stick

- Easily identify and remove all infected dentin.
- Single-patient applicators make tray set-up easy.
- No syringe, no sheath, no tips to replace.
- No cross-contamination, no messy spills, no stains!
- Made in the USA.



IMPRESSIONING

Versatile choices for easier, safer impressions

UNIQUE AND COST-EFFECTIVE SOLUTIONS FOR ALL IMPRESSION TECHNIQUES



NoCord™ VPS

One-Step, Self-Retracting Impressioning System

- The safest, most unique VPS available today.
- Built-in hemostasis – wash stops bleeding and fluid flow for clear margins, and provides exceptional detail.
- MegaBody™ Tray Material creates sulcular expansion. Enables complete margin capture, usually without retraction cord. Adds extra hemostasis when you do use cord.



Made in the USA

Access GOLD®

- Low-viscosity wash material, syringes easily into sulcus and around prep.
- High-contrast color makes it easy to verify all margins have been captured.

Access MONOPHASE®

- Non-slumping, medium-body material.
- Perfect for edentulous impressions or as a tray material.



Access HEAVY®

- High-viscosity tray material, and the perfect complement to Access GOLD.
- Easy to place, easy to remove.

Access PUTTY®

- Non-oily and easy to mix.
- For use as a tray material or for easy fabrication of stents, matrices and more.

EASIER, FASTER BITE REGISTRATION



Access BLUE®

VPS Quick Bite Registration Material

- Mousse-like viscosity for easy, fast application.
- Sets hard, like laboratory plaster in just 30 seconds.
- Most accurate bite registration for intercuspal and centric relationships.
- Cut and trim the impression without it crumbling.
- Made in the USA.



3.n.One™

Impression Trays For Dual-Arch Technique

- Unique undercut ridges lock impression materials in place.
- High rigidity with non-tear, flexible webbing ensures maximum accuracy.
- Three different shapes to meet all applications.

BETTER RETRACTION WHEN YOU NEED IT



Access® Edge

High-Viscosity Kaolin Clay with 15% Aluminum Chloride

- A thicker, more viscous paste that provides more 'push' of the soft tissue.
- Uses standard composite syringe; no glass vials.
- Dual action – stops bleeding, and absorbs moisture.
- Made in the USA.
- GingiCap™ Compression Caps included.

Two Kaolin Clay formulations to meet every need!



AccessFLO®

Low Viscosity Kaolin Clay with 15% Aluminum Chloride

- Flowable paste that more easily penetrates the sulcus.
- Dual action – 15% aluminum chloride plus kaolin clay stops bleeding and absorbs moisture.
- Made in the USA.
- GingiCap™ Compression Caps included.



GingiTrac™

VPS Gingival Retraction System

- Makes digital dentistry easier and better.
- Unique VPS material gently retracts tissue yet removes cleanly with no debris or particles that can distort digital images.
- Ammonium aluminum sulfate to control bleeding and oozing.
- Made in the USA.
- GingiCap™ Compression Caps included.

POST-ENDO TREATMENT

Centrix products provide a superior approach

ROOT CANAL TREATMENT



Calciject®

Calcium Hydroxide Paste for Root Canal Therapy

- No messy mixing of powder and liquid.
- Direct injection from package to mouth.
- Single-dose tips for maximum asepsis.
- Bendable 20 ga NeedleTube® tips for easy canal access.
- Superior results compared to spiral paste fillers.
- Bacteriostatic and radiopaque.
- Made in the USA.



Connexio®

Dual-Cure, Self-Etch Bond for Core Buildups

- Designed for the demands of deep core build-ups where light cannot reach.
- Dual-cure adhesive system provides the highest possible bond strength.*
- Nanoparticle reinforced.
- Water solvent (not acetone or ethyl alcohol); makes Connexio more tolerant of humid canal conditions than other adhesives.

*Migliau, G; Piccoli, L; Di Carto, et al, Comparison between three glass fiber post cementation techniques. (Ann. Stomatol (Roma); 2017 Jan-Mar: 8(1): 29-33. Published online 2017 Jul 3 doi [10.11138 ads/2017.8.1.029].

DUAL-CURE CORE BUILDUPS THAT BOND

More dentists choose **Encore D/C** and **Connexio** for a bond that can't be broken!



Encore® D/C MiniMix®

Dual-Cure, Core Buildup Composite

- Non-slumping formula stays where it's placed, will not run and does not require a matrix or core form.
- Exclusive MiniMix® cartridge (0.5 mL) controls costs while eliminating waste and cross-contamination.
- Cuts like dentin to reduce ditching.
- NeedleTube® mixing tip delivers material to bottom of the post hole.
- Made in the USA.
- Also available in traditional 25 mL automix cartridges.



ESSENTIAL TOOLS FOR STRUCTURAL RESTORATIVE INTEGRITY




AbsoLute™

Universal Resin Cement

- Self-etch, self-adhesive, dual-cure for faster, easier placement. No primer or conditioner needed.
- Definitive cementation that lasts.
- Designed especially for feldspathic, high leucite, lithium disilicate, alumina and zirconia materials.
- Made in the USA.

3 Shades = A natural-looking restoration

- **Translucent:** Provides a natural esthetic
- **Opaque (White):** Completely blocks underlying colors/show-through
- **Universal (Opaque Tooth Shade):** Completely blocks underlying colors/show-through



Onyx™

40% Phosphoric Acid Etching Gel

- Highest strength, most effective etchant.
- Black silicon carbide particles provide superior contrast.
- Better control and easier placement.
- Stays where placed, will not run.
- Particles can be agitated on tooth to mechanically scrub and clean the surface.
- No dye, does not stain, rinses completely.
- Made in the USA.

TEMPORIZATION

Complete solutions designed to work better together

THE RIGHT FILLING MATERIAL FOR EVERY CLINICAL SITUATION



Short-Term
Temps

Tempit® Moisture-Activated Temporary Filling & Sealing Material

- Ideal for endo access openings and short-term temporaries.
- Moisture-activated, sets in 5 minutes or less.
- Expands slightly when set, seals out microleakage and bacterial ingress.
- Contains calcium sulfate and zinc oxide.
- Significant bacteriocidal activity up to 14 days.¹

¹ Antibacterial Properties of Temporary Filling Materials, JOE, Volume 32, Number 3, March 2006



Short-Term Temps
with Eugenol

Tempit®-E Moisture-Activated Temporary Filling & Sealing Material with Eugenol

- Palliative eugenol formulation for temporary filling use, especially in cases of reversible pulpal inflammation.
- Anti-inflammatory and local anesthetic effects helps facilitate pulpal healing.
- Moisture activated, sets in five minutes.
- Expands slightly when set to seal out microleakage and bacterial ingress.
- Can be used as a liner in deep preparations.



Indirect Temp
Restorations

Tempit® L/C Light-Activated Temporary Material for Inlays & Onlays

- Perfect for indirect restorations, inlays/onlays.
- Seals implant screw access openings.
- Sets firm yet flexible.
- Easily removed with an explorer without damaging the margins.
- Microfilled polymer resin.
- Light-activated, sets in 40 seconds or less.
- No etching or bonding agent required.



Longer-Term
Temps

Tempit® Ultra-F Light-Activated, Longer-Term Temporary Material

- Perfect for longer term temporaries.
- Rigid, harder wearing material, with fluoride release.
- High compressive strength for bruxism patients.
- Ideal for tooth-colored restorations in deciduous teeth.
- Less soluble than reinforced ZOE intermediate restorative material.
- Light-activated, sets in 40 seconds or less.
- No etching or bonding agent required.
- **Tempit Ultra: without fluoride**



THE DENTAL
ADVISOR
++++½

TWO BETTER TEMPORARY CEMENTS

THE DENTAL
ADVISOR
++++½



Actual size.

Access® Automix Temporary Cement

- Perfect for longer-term temporary use.
- Low film thickness.
- Flows and places easily.
- Semi-elastic set.
- Safe, easy adhesive.
- Calcium hydroxide provides bacteriostatic effect on mold.
- Eugenol-free, won't interfere with bonding agents or other resin materials.

NoMIX® Short-term, Temporary Cement

- Easy to use, dispense, clean and remove from prep.
- Low film thickness, harder set.
- Tube applicators are ideal for both chair-side application and for patient take-home.
- Expands slightly upon setting for a tight margin seal.
- Moisture activated, no need to dry.
- **Also available in syringe delivery for chair-side application.**



THE DENTAL
ADVISOR
++++½




Access Crown® Temporary Crown & Bridge Material

- **More:** shades (7), speed, strength and savings.
- **Ideal handling characteristics:** ready to remove from mouth in 60 seconds.
- Low shrinkage, low heat generation.
- High compressive, flexural strength.
- Highly aesthetic – polishes easily, beautifully.




Natural Esthetics, High Fluorescence, Polishes Easily.

✓ MORE SPEED Intra Oral Set & Final Cure Time



✓ MORE STRENGTH



DELIVERY SYSTEMS

Control costs and eliminate cross-contamination

APPLY MATERIALS EASILY AND MORE PRECISELY

Benda® Brush

The Original Precision Disposable Brush Applicator

- When you need to apply more than a drop.
- Bend head to any angle with no bounce-back!
- Unique non absorbing fibers have a greater carrying capacity, ensuring better spreading.
- Regular point for large applications.
- Mini point for smaller spaces.
- Made in the USA.



Benda® Micro

Precision Disposable Fiber Applicator

- Perfect for when you need only a drop.
- Regular point (2.0 mm).
- Fine point (1.5 mm).
- Extra-Fine point (1.0 mm).
- Non-linting design will not leave fibers behind in restorations or sulcuses.
- Made in the USA.

Benda® Twin

Double-ended, Twin-style Applicator

- When you are using more than one material.
- Ideal for 2-step procedures.
- Extended reach.
- Made in the USA.

Benda Twin with snap-apart handle. The lowest cost alternative!

Benda Brush

Benda Micro

More sizes, more styles and more colors means the right applicator for every material and every need!



DELIVER MATERIALS MORE EFFECTIVELY



More Usable Material Per Cartridge!

■ Mixer Waste (g)
■ Usable Material (g)



Less waste in every tip!

SuperMixer™

AutoMix Mixing Tips

- High-efficiency mixer, up to 38% reduced material waste.
- Up to 50% more units from every cartridge.
- More consistent mixing of any material.
- Color-coded hub for easy replacement of old-style tips.
- Replaces any 1:1 mixing tip – same color for same color.

AccuDose® & NeedleTube®

Delivery Systems

- Dispense exactly the material needed.
- Prevent cross-contamination.
- Ready to use, dispense directly into the prep.
- No waste, no contamination.
- No more 'composite keeper'
- Compatible with all C-R Syringes and guns.
- Made in the USA.

PREVENTION

FLUORODOSE® UNIT DOSE, 120-CT.

- 360086 Mint, 0.3 mL LoliTray™
- 360087 Bubble Gum, 0.3 mL LoliTray™
- 360107 Melon, 0.3 mL LoliTray™
- 360110 Cherry, 0.3 mL LoliTray™
- 360137 Caramel, 0.3 mL LoliTray™

CHAMP®

- Unit Dose**
- 360116 16-ct., 0.25 mL LoliPack
- 360120 64-ct., 0.25 mL LoliPack
- 310033 Champ Sealant Syringe Kit 4x1.2 mL, 20 tips
- 310034 Champ Sealant Syringe Single 1x1.2 mL, 5 tips
- 310032 22 ga Luer Lock Needle Tips, 20-ct.
- 640001 22 ga Luer Lock Needle Tips, 100-ct.

LOLLICAINE® UNIT DOSE, 120-CT.

- 360088 Bubble Gum, 0.3 mL LoliPack
- 360090 Cherry, 0.3 mL LoliPack
- 360092 Mint, 0.3 mL LoliPack
- 360094 Piña Colada, 0.3 mL LoliPack

E-LIEVE™

- 310171 50 x 0.5 mL Applicators

EXPOSÉ®

- 310334 Fine, 100-ct., Blue
- 310330 Fine, 576-ct., Blue
- 310335 Extra Fine, 100-ct., Blue
- 310331 Extra Fine, 576-ct., Blue

IMPRESSIONING

NOCORD™ VPS

- 310280 VPS Intro Trial Kit
- 310285 VPS System Intro Kit
- 310282 Wash Material, 2-ct.
- 310284 Wash Material, 10-ct.
- 310281 MegaBody™ Tray Material, 2-ct.
- 310283 MegaBody Tray Material, 10-ct.
- 310286 MegaBody, 380 mL, 1-ct.
- 310287 MegaBody, 380 mL, 2-ct.

ACCESS®

- 360002 GOLD®, 2 x 50 mL
- 360016 HEAVY®, 2 x 50 mL
- 360017 MONOPHASE®, 2 x 50 mL
- 360015 PUTTY®, 2 x 280 mL Jars

ACCESS RETRACTION PASTE

- 360113 Access® Edge, 15-ct.
- 360114 Access® Edge, 60-ct.
- 640002 16 ga Luer Lock Needle Tip, 30-ct.
- 360122 AccessFLO®, 15-ct.
- 360123 AccessFLO®, 60-ct.

GINGITRAC®

- 360100 Refill Kit, 2 x 50 mL
- 360101 Standard Kit, 2 x 50 mL, 1:1/2:1 Gun
- 360102 Plus Kit, 4 x 50 mL
- 360068 MiniMix Starter Kit, 24-ct.
- 360069 MiniMix Refill Kit, 48-ct.

ACCESS® BLUE

- 360006 Bite Registration Kit, 2 x 50 mL

3.N.ONE™ IMPRESSION TRAYS

- 360011 Anterior, 35-ct.
- 360012 Posterior, 50-ct.
- 360013 Quadrant, 35-ct.
- 360014 Assorted, 40-ct.

POST-ENDO

CALCIJECT®

- 320003 Prefilled Tip, 0.20 g, 24-ct.

CONNEXIO®

- 310270 Self-Etch Bond, 3 mL Part A & 3 mL Part B

ENCORE® D/C

MiniMix® Cartridge 0.5 mL

- 310129 Natural, 12-ct.
- 310254 Natural, 36-ct.
- 310130 Contrast Blue, 12-ct.
- 310255 Contrast Blue, 36-ct.
- 700000 MiniMix Kit, 12-ct., Natural, with Snub-Nose Snap-Fit Syringe

Automix Cartridge 25 mL, 10 mix tips

- 310150 Natural, 1 x 25 mL
- 310152 Contrast Blue, 1 x 25 mL

MiniMix® Nozzles

- 910630 AutoMix 1:1 Mixing Nozzles, 10-ct.

ABSOLUTE® UNIVERSAL RESIN CEMENT

- 310113 Opaque, 1 x 5 g Syringe
- 310114 Translucent, 1 x 5 g Syringe
- 310015 Universal, 1 x 5 g Syringe
- 310016 Assorted, 3 x 5 g Syringe
- 210008 Blunt Mixing Tips, 50-ct.
- 210007 Tapered Mixing Tips, 50-ct.
- 210009 IntraOral Tips, 50-ct.

ONYX® ETCH

- 310053 6 x 1.2 mL Syringes, 30 Needle Tips
- 310054 15 x 1.2 mL Syringe, 50 Needle Tips
- 310055 15 x 1.2 mL Syringe, 100 Needle Tips
- 310052 25 x 1.2 mL Syringe Only
- 310032 22 ga Luer Lock Needle Tips, 20-ct.
- 640001 22 ga Luer Lock Needle Tip, 100-ct.

TEMPORIZATION

TEMPIT® LINE, PREFILLED TIPS

- 310060 Tempit, 0.35 g, 30-ct.
- 310031 Tempit-E, 0.35 g, 30-ct.
- 310063 Tempit L/C, 0.25 g, 30-ct.
- 310064 Tempit Ultra-F, with fluoride, 0.20 g, 30-ct.
- 310065 Tempit Ultra, no fluoride, 0.20 g, 30-ct.
- 310067 Tempit Essentials Kit (5 ea. of Tempit, Tempit-E, Tempit L/C, Tempit Ultra-F)

ACCESS® CROWN®, 76 g CARTRIDGE

- 360029 Shade A3.5
- 360030 Shade A3
- 360031 Shade B1
- 360034 Shade A2
- 360035 Shade A1
- 360037 Shade C2
- 360038 Shade BL

ACCESS® TEMPORARY CEMENT

- 310521 1 x 5 mL Automix Cartridge, Mixing Nozzles, 10-ct.

NOMIX® TEMPORARY CEMENT

- 310520 Starter Kit, 1 mL prefilled syringe, 6-ct., Take-Home Applicator, .50 g, 24-ct., 30 Syringe Tips
- 310157 Take-Home Applicator, .50 g, 50-ct.
- 290040 Syringe Refill Kit, 12 x 1 mL prefilled syringe, 60 Syringe Tips
- 640002 16 ga Luer Lock Needle Tip, 30-ct.

DELIVERY SYSTEMS

ACCUDOSE® & NEEDLETUBE®

- 290032 19 ga NeedleTube, 100-ct.
- 220008 White Plug, 100-ct.

- 290031 20 ga NeedleTube, 100-ct.
- 220000 Green Plug, 100-ct.

- 290049 24 ga NeedleTube, 100-ct.

- 920201 AccuDose TriPack Starter Kit, 75 Tubes & Plugs

- 210021 Low-Viscosity Tube, 100-ct.
- 210023 Low-Viscosity Tube, 500-ct.
- 210000 Gray Plug, 100-ct.
- 280021 High-Viscosity Tube, 100-ct.
- 280023 High-Viscosity Tube, 500-ct.
- 290046 Black Plug, 100-ct.

- 210001 Separate Clear Tube, 100-ct.
- 210003 Separate Clear Tube, 500-ct.
- 240001 Separate Black Tube, 100-ct.
- 250011 PCR Clear Tube, 100-ct.
- 250001 PCR Opaque Tube, 100-ct.
- 310032 22 ga Luer Lock Needle Tip, 20-ct.
- 640001 22 ga Luer Lock Needle Tip, 100-ct. for Onyx & Champ Syringes
- 640002 16 ga Luer Lock Needle Tip, 30-ct. for Access Edge & NoMIX
- 290038 22 ga Needle Tip, Gray, 144-ct.

SUPERMIXER™ MIXING NOZZLES

- 910628 Yellow hub, 50-ct.
- 910631 Yellow hub, 18 ga tip, 50-ct.
- 910632 Pink hub, 50-ct.
- 360010 Teal hub, 24-ct.
- 910626 Blue hub, 50-ct.
- 910630 AutoMix 1:1 Mixing Nozzles, 10-ct.

BENDA® BRUSH – BRUSH APPLICATORS

- 370101 Regular, Assorted, 144-ct.
- 370501 Regular, Red, 576-ct.
- 370502 Regular, Blue, 576-ct.
- 370503 Regular, Yellow, 576-ct.
- 370504 Regular, Green, 576-ct.

- 370111 Mini Assorted, 144-ct.
- 370601 Mini, Red, 576-ct.
- 370602 Mini, Blue, 576-ct.
- 370603 Mini, Yellow, 576-ct.
- 370604 Mini, Green, 576-ct.

- 370105 Twin, Pink, 200-ct.
- 370106 Twin, Blue, 200-ct.
- 370107 Twin, Pink, 400-ct.
- 370108 Twin, Blue, 400-ct.

BENDA® MICRO – FIBER APPLICATORS

- 380005 Micro Regular, Assorted, 400-ct.
- 380014 Micro Regular, Blue, 576-ct.
- 380015 Micro Regular, Green, 576-ct.
- 380016 Micro Regular, Purple, 576-ct.
- 380017 Micro Regular, Silver, 576-ct.

- 380006 Micro Fine, Assorted, 400-ct.
- 380012 Micro Fine, Yellow, 576-ct.
- 380013 Micro Fine, Red, 576-ct.

- 380033 Micro XF, Assorted, 400-ct.
- 380029 Micro XF, Light Green, 576-ct.
- 380030 Micro XF, Ivory, 576-ct.

- 380021 Micro Twin Regular, Blue, 200-ct.
- 380023 Micro Twin Regular, Blue, 400-ct.
- 380020 Micro Twin Regular, Pink, 200-ct.
- 380022 Micro Twin Regular, Pink, 400-ct.
- 380028 Micro Twin XF, Green, 400-ct.
- 380027 Micro Twin XF, Ivory, 400-ct.



FREE LUNCH & LEARN

- Easy. Esthetic. Durable. Temporization.
- Prevention for Life®
- Impressioning
- The Core Issue

To schedule a Centrix Lunch & Learn, call 800.235.5862 or visit centrixdentallearning.com

CENTRIX INCORPORATED

770 River Road | Shelton, CT | 06484 | USA

203.929.5582, Toll-Free (USA & Canada): 800.235.5862

Centrix is an ISO 13485 certified dental manufacturer and our products are marketed world-wide with FDA, CE, Health Canada, TGA and PAL registrations on file. © 2021 Centrix, Inc. | #GUIDES, REV D



centrix®
MAKING DENTISTRY EASIER™